orthogum

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

Polysiloxane, condensation-reaction silicone-elastomer. Type 0 mixture, (ISO 4823 / ADA 19: **Type 0 - Putty Consistency**).

Product identification: Particularly soft and elastic mass both before and after polymerisation. Indicated for orthodontics and complete dentures.

Trays: all rigid types; in particular stainless steel trays (Hi-Tray Zhermack); for trays without holes and retention trays, use the special adhesive (Universal Tray Adhesive Zhermack).

Warning: Make sure the catalyst tube or indurent liquid bottle are securely closed after use. The silicone-elastomers are chemically resistant, therefore avoid staining overalls or clothers. Use gloves.


Warning: In sensitive subjects, polysiloxanes may cause irritation or other allergic reactions.

Cleaning and disinfection of the impression: After rinsing well under running water, the impression can be immediately disinfected. Ideal disinfecting is achieved by immersing the impression in **Zeta 7 solution** (Zhermack) or spraying it with **Zeta 7 spray** (Zhermack). When using other disinfectants, always follow the manufacturer's instructions.

Storage: Orthogum impressions should be stored at room temperature. Orthogum impression material is guaranteed for a period of 24 months from the date of manufacturing, if stored correctly at 5° to 27°C (41°-80° F).

Models: The impression should be poured in a period of from 30 minutes to 72 hours after polymerization. In order to avoid surface tension use Zhermack TENSILAB. Orthogum is compatible with all high quality plaster products.

Electroplating: Orthogum impressions can be silver or copper plated.

*The above mentioned times are intended from the start of mixing at 23°C - 73°F.

Important remarks: Consulting services on the use of our products, provided in any form, be that verbal, written or through demonstrations, are based upon the current state of dentistry and of our know-how. They are to be considered as non-binding information, also in relation to possible third party claims, and do not exonerate the user from personally checking whether the product is suitable for the planned use. Use and application by the user take place under no possible control by the Company, and are therefore under the user's responsibility. Possible liability for damages is limited to the value of goods supplied by the Company and utilized by the user.
FOR DENTAL USE ONLY

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