

Pulpdent Calcium Hydroxide Pastes

Multi-Cal™

MULTI Multi-Cal Kit, 4 x 1.2 mL syringes
+ 8 applicator tips (22 gauge x ½")

MULTI-3 Multi-Cal 3 mL syringe



40% ($\pm 2\%$) calcium hydroxide in a creamy aqueous ethylcellulose paste

Indicated for root canal therapy and vital pulp therapy

Multi-Cal is a non-setting, pre-mixed paste easily removed from root canals with file and irrigation. Radiopaque.

Packaged in 1.2 mL or 3 mL push syringes for direct dispensing through 22-gauge x ½" applicator tips.

World Renowned

Since developing the first pre-mixed calcium hydroxide-aqueous methylcellulose pulpal dressing in 1947, the Pulpdent name has been synonymous with calcium hydroxide.

Multi-Cal, TempCanal, Pulpdent Paste are water based pastes that dry with air but do not set hard. Forendo Paste is silicone oil based and neither dries nor sets hard. All four pastes are easily removed from the canal with file and irrigation.

- **Non-setting, pre-mixed pastes**
- **Bactericidal and bacteriostatic**
- **Radiopaque**
- **pH > 12**

Forendo Paste™

FORE Forendo Paste, 2.2 gm syringe
+ 20 applicator tips



22% calcium hydroxide with iodoform in a silicone oil-based paste

Indicated for root canal therapy.

Iodoform combined with calcium hydroxide provides a strong antibacterial dressing for root canal therapy. Packaged in a 2.2 gram push syringe with applicator tips for direct dispensing into canals. A non-setting, pre-mixed paste easily removed from root canals with file and irrigation. Radiopaque.